## Linola® sept Mund- und Rachenspülung Linola® sept Mouth and Throat Rinse

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Januar 2021

#### Dear Patient,

**Please read the following package leaflet carefully,** it contains important information that you need to be aware of when using this medical device. Please keep this leaflet. You may wish to read it again later. If you have any further questions, consult your doctor or pharmacist.

### Linola® sept Mouth and Throat Rinse UK – Instructions for Use –

**Corona prophylaxis by physical reduction of viral load in the mouth and throat.**

#### Indications for use

Linola sept Mouth and Throat Rinse

- for physical cleansing
- for reducing the viral load in the oral cavity and throat.

#### Product information

Coronaviruses and other pathogens penetrate the body via the mouth and throat and reproduce quickly. They can then spread in the body and be transmitted to other people through droplets when talking, sneezing or coughing.

The surfactants contained in Linola sept Mouth and Throat Rinse adhere to the outside of the viruses. This means that the viruses cannot bind to the human cells and thus cannot infect the mouth and throat.

By gargling Linola sept Mouth and Throat Rinse, the viral load in the oral cavity and throat is reduced, minimising the risk of coronavirus transmission. The product is therefore intended for corona prophylaxis.

Linola sept Mouth and Throat Rinse does not provide protection against infection with coronaviruses on its own. We recommend using it to help with corona prophylaxis in addition to the well-known measures (personal hand hygiene, wearing a mask over the mouth and nose, social distancing, frequently ventilating rooms). Clinically tested.

#### Application

Shake Linola sept Mouth and Throat Rinse gently and use the dosing cap to measure approximately 10–20 ml.

Shaking mixes the ingredients in the mouth rinse causing the product to foam up.

Do not dilute the rinse.

Gargle in the throat and swirl around the mouth for 60 seconds.

Do not swallow the product.

Do not rinse out the mouth and throat with water.

Clean the dosing cap after use.

Use Linola sept Mouth and Throat Rinse 2 to 3 times a day, especially before and after contact, to reduce the viral load.

#### Possible side effects

Temporary symptoms such as tingling or burning on the tongue or in the mouth can occur in rare cases.

If you notice any other side effects, please consult your doctor or pharmacist.

#### Restrictions for use

In case of allergies or suspected allergies to any of the ingredients, Linola sept Mouth and Throat Rinse must not be used.

#### Important information

Can be also used by children over 6 years.

Does not contain alcohol.

The mouth and throat rinse must not be used in children and individuals who may not have control over their swallowing reflex.

#### Storage and shelf life

Store Linola sept Mouth and Throat Rinse out of the sight and reach of children.

Do not store above 25°C.

Shelf life once opened: 12 months.

The expiry date is printed on the flap of the outer packaging and/or on the label.

The expiry date refers to the last day of that month.

Do not use Linola sept Mouth and Throat Rinse after this date.

#### Pack sizes

Linola sept Mouth and Throat Rinse is available as a 375 ml pack.

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Dr. August Wolff GmbH & Co. KG Arzneimittel